Quit DNAnuclospin Plant II

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100

1

Table 1: Kit specifications at a glance

Parameter	NucleoSpin® Plant II
Technology	Silica-membrane technology
Format	Mini spin columns
Sample material	Up to 100 mg wet weight Up to 20 mg dry weight
Lysate clarification	NucleoSpin® Filters
Fragment size	50 bp– approx. 50 kbp
Typical yield	1–30 µg
A_{260}/A_{280}	1.8–1.9
Elution volume	2 x 50 µL
Preparation time	30 min/prep
Binding capacity	50 µg

Table 2: Plant species tested with NucleoSpin® Plant II

Plant species	Plant tissue / organ	Lysis buffer successfully tested	
		PL1	PL2
<i>Abies alba</i> (fir)	Needle	✓	✓
<i>Acer griseum</i>	Leaf	✓	✓
<i>Amorphophallus titanum</i>	Leaf	✓	Not tested
<i>Apium graveolens</i> (celery)	Corm	✓	✓
<i>Arabidopsis thaliana</i>	Leaf	✓	✓
<i>Boreava orientalis</i>	Leaf, herbarium sample	✓	✓
<i>Carex annectens</i>	Leaf	✓	✓
<i>Carex waponahkikensis</i>	Leaf, silica-gel dried	✓	✓
<i>Cleisostoma racemiferum</i>	Inflorescence rachis, silica-gel dried	✓	Not tested
<i>Doritis pulcherrima</i>	Leaf, silica-gel dried	✓	Not tested
<i>Eichornia azurea</i>	Leaf	✓	Not tested
<i>Encephalartos natalensis</i>	Leaf	✓	Not tested
<i>Galium aparine</i>	Leaf	✓	✓
<i>Hordeum</i> sp. (barley)	Leaf	✓	✓
<i>Isatis kotschyana</i>	Leaf, herbarium sample	✓	✓
<i>Laurus azorica</i> (laurel)	Leaf	✓	Not tested
<i>Lupinus</i> sp. (lupin)	Leaf	✓	✓
<i>Lycopersicon esculentum</i> (tomato)	Stem	✓	✓
<i>Myagrum perfoliatum</i>	Leaf, herbarium sample	✓	✓
<i>Oryza sativa</i> (rice)	Leaf	✓	✓
<i>Persea feru./caerulea</i>	Leaf	✓	Not tested
<i>Pteridium</i> sp.	Leaf	✓	Not tested
<i>Pterocarya fraxinifolia</i>	Leaf	✓	Not tested
<i>Quercus cerris</i>	Leaf	✓	✓
<i>Quercus frainetto</i>	Leaf	✓	✓
<i>Rosa</i> sp. (rose)	Leaf	✓	✓
<i>Rubus fruticosus</i> (blackberry)	Leaf	✓	✓
<i>Sameraria nummularia</i>	Leaf, herbarium sample	✓	✓
<i>Secale</i> sp. (rye)	Leaf	✓	✓
<i>Stereochilus</i> sp.	Leaf, silica-gel dried	✓	Not tested
<i>Tauscheria lasiocarpum</i>	Leaf, herbarium sample	✓	✓
<i>Trachycarpus takil</i>	Leaf	✓	Not tested
<i>Trichoglottis</i> sp.	Leaf, silica-gel dried	✓	Not tested
<i>Triticum aestivum</i> (wheat)	Leaf	✓	✓
<i>Vigna radiata</i> (mung bean)	Root	✓	✓
<i>Zea mays</i> (maize)	Leaf	✓	✓
<i>Zea mays</i> (maize)	Grain, dried, ground coarsely	✓	✓
Fungal mycel (not specified)		✓	Not tested

NucleoSpin® Plant II**50 preps****REF 740770.50**

Wash Buffer PW2
(Concentrate) 25 mL
add 100 mL
ethanolRNase A 6 mg
dissolve in
600 µL H₂O